

Exp No :

Date:

Refractive Index of the Material of a Glass slab

Aim: To determine the refractive index of a glass slab using a Traveling Microscope.

Apparatus: Glass Slab, Traveling Microscope, Lycopodium powder, paper etc.

Theory:

$$\text{Refractive index of the material of the glass slab} = \frac{\text{Real Thickness of the Glass Slab}}{\text{Apparent Thickness of the Glass Slab}}$$

Observations:

Value of One main scale division 1 MSD = mm

Number of Vernier Scale Divisions N = div

Least Count LC = **1 MSD – 1 VSD**

$$= \frac{1 \text{ MSD}}{N} = \text{mm}$$

$$= \text{cm}$$

$$\text{Total Reading} = \text{MSR} + (\text{VSR} \times \text{LC})$$

Reading	Reading of the vertical scale of the Microscope when focused on the			Real Thickness $R_3 - R_1$	Apparent Thickness $R_3 - R_2$	n = $\frac{R_3 - R_1}{R_3 - R_2}$
	Cross-mark without Glass Slab	Cross-mark with Glass Slab	Lycopodium Powder on Glass Slab			
MSR (cm)						
VSR						
Total (cm)	R ₁ =	R ₂ =	R ₃ =			

$$\text{Refractive index of the material of the glass slab} = \frac{\text{Real Thickness}}{\text{Apparent Thickness}}$$

$$n = \frac{R_3 - R_1}{R_3 - R_2} =$$

Result:

$$\text{Refractive index of the material of the glass slab} =$$